

ReadSense AI — Pipeline Architecture Explained

What this is: A plain-language breakdown of exactly what ReadSense does, why we use Whisper + Librosa + Gemma together, and how data flows from microphone to diagnosis. Written for the team.
Last updated: 2026-07-14

What ReadSense Is Actually Trying to Do

ReadSense is a **reading diagnostics tool for children**. When a child reads a sentence out loud, a teacher cannot listen to 30 kids at once and individually diagnose *why* each one is struggling. ReadSense solves this.

It does **three distinct things**:

1. Measure *How* They Read (Acoustic Analysis)

Not just *what* they said — but *how* they said it:

- **How fast?** (Words Per Minute)
- **Where did they pause?** (and for how long?)
- **Did their voice pitch spike?** (anxiety/frustration signal)
- **Did they repeat words or stumble?** (hesitation markers)

2. Diagnose *Why* They Struggle (5 Gemma Agents)

Five specialized AI agents receive that acoustic data and reason about it:

Agent	Job
Phonetic Diagnostician	Which specific sounds/phonemes are hard? ("br" blends, "th" sounds)
Difficulty Classifier	Is the struggle phonetic, cognitive, fluency-based, or mixed?
Engagement Monitor	Is the child anxious, frustrated, or confident right now?
Practice Generator	Generate 5 personalized practice sentences targeting their weak spots
Progress Tracker	Compare today to past sessions — are they improving?

3. Report to the Teacher

A clean dashboard showing: score, struggling words, emotional state, and what to practice next. Built for teachers who have no time to manually analyze each child.

What the Libraries Actually Do

These are **not** just "audio to text" tools. They extract precise diagnostic measurements.

[openai-whisper](#) — Word-Level Transcription with Timestamps

Whisper returns **timestamps for every single word**, not just the full transcript:

```
result = whisper.transcribe("reading.wav")

# What you actually get back:
{
  "text": "The brave dog ran quickly",
  "segments": [
    {"word": "The", "start": 0.00, "end": 0.18, "probability": 0.99},
    {"word": "brave", "start": 0.18, "end": 0.52, "probability": 0.97},
    {"word": "dog", "start": 0.52, "end": 0.73, "probability": 0.99},
    {"word": "ran", "start": 1.20, "end": 1.45, "probability": 0.95}, # ← gap
    before 'ran' = pause!
    {"word": "quickly", "start": 1.45, "end": 1.90, "probability": 0.91},
  ]
}
```

From this alone you can calculate:

- **WPM** — $(5 \text{ words} \div 1.90 \text{ seconds}) \times 60 = 158 \text{ WPM}$
- **Exact pause locations** — the gap between "dog" ending at 0.73s and "ran" starting at 1.20s = **470ms pause**
- **Word confidence** — low **probability** means the child mumbled or mispronounced that specific word

librosa — Acoustic Signal Analysis


Whisper gives you the *linguistic* picture. Librosa gives you the *physical* sound picture:

```
What librosa measures:
├─ Pitch (F0)          - Is the child's voice rising? Falling? (anxiety = rising pitch)
├─ Pitch Variance      - How much does pitch fluctuate? (high variance = stress/nervousness)
├─ Silence Intervals  - Detect pauses whisper might miss (sub-200ms micro-hesitations)
└─ Energy/Volume       - Is the child reading louder at certain words? (effort/struggle signal)
```

The Full Pipeline — Step by Step

Here is exactly what happens when a child reads a sentence:

```
Child speaks into mic
  ↓
Browser captures audio (WebM blob via MediaRecorder API)
  ↓
Frontend sends audio_base64 to FastAPI backend
```



✓ Why This Approach Is Better Than Feeding Raw Audio to Gemma

We tested feeding raw audio directly to `gemma-4-31b-it` and it returned:

```
400 Bad Request: Audio input modality is not enabled for this model
```

This is because the **31B Gemma model does not support audio**. Only the smaller E2B, E4B, and 12B variants do — and those are not available via the Google AI Studio API, only as local downloads.

The Doctor Analogy

Imagine you're a specialist doctor (the Gemma agent). Would you rather:

- **Option A:** Receive a chart saying: "*Patient heart rate: 142bpm, blood pressure: 150/95, temperature: 38.5°C*" → You immediately reason and diagnose.
- **Option B:** Receive a raw EKG printout and be asked to extract all the numbers yourself before diagnosing.

Whisper + Librosa are your lab instruments. They produce precise measurements. **Gemma is the specialist doctor** who receives clean, structured data and makes an expert judgment.

The 31B Gemma model is optimized for **reasoning and language understanding**, not raw audio signal processing. Giving it pre-extracted acoustic features makes its 5 diagnoses **more accurate and more reliable** than raw audio would.

Gemma Model Clarification (Important)

The blueprint has a typo (`gemma-4-9b-it` — this model does not exist). Here is the correct model lineup:

Model	Audio via API	Best For
<code>gemma-4-e2b-it</code>	✗ Local only	Mobile/Edge
<code>gemma-4-e4b-it</code>	✗ Local only	Mobile/Edge
<code>gemma-4-12b-it</code>	✗ Local only	Audio (but no API)
<code>gemma-4-26b-it</code> (MoE)	✓ AI Studio	Text reasoning
<code>gemma-4-31b-it</code> (Dense)	✓ AI Studio	Text reasoning — this is our model

Our architecture is correct: use `gemma-4-31b-it` via the API for all 5 reasoning agents, and use Whisper+Librosa locally for audio preprocessing.

Key Takeaway

- ReadSense does not ask Gemma to *listen* to a child read.
It asks Gemma to *diagnose* a child's reading based on precise acoustic measurements.
Whisper and Librosa are the instruments. Gemma is the expert that interprets the results.

This separation is what makes the pipeline fast, reliable, and accurate — and it's exactly what the original blueprint intended.